

TiNoI-1TM

NITINOL SPRINGS

PREMIUM NITINOL

by Modern Orthodontics



Open or Closed
Constant Force
Diamond Polished
Medium Force





TiNol-1TM Nitinol Springs

TiNol-1TM NiTi open coil spring: Open coil springs **open and maintain spaces with constant, gentle force**. From start to finish this spring delivers constant, unvarying force for predictable results-even when spring compression is changed the force stays the same.

.010X.030 spring in 15" lengths

TiNol-1TM NiTi open coil spring

.010" x .030" 15" on spool NT1030O-15I

.012" x .030" 15" on spool NT1230O-15I

TiNol-1TM NiTi closed coil spring with eyelets:
Use closed coil eyelet springs to close spaces with ease, **without patient discomfort** and within a predictable time period.

.010X.030 spring in 9mm and 12mm lengths
in packs of 10

TiNol-1TM NiTi closed coil spring

.010" x .030" 9mm pk of 10 NT1030C-09M

.010" x .030" 12mm pk of 10 NT1030C-12M

TiNol-1TM NiTi closed coil spring with eyelets for class II correction:
Use it in place of class II elastics when your patient compliance is a problem. Tie it once and it would work permanently for months. Spring should be covered with soft sleeve to eliminate pinching and food trapping.

TiNol-1TM .010X.030 spring in 15 mm and 18 mm lengths
in packs of 10

TiNol-1TM NiTi closed coil spring

.010" x .030" 15mm pk of 10 NT1030C-15M

.010" x .030" 18mm pk of 10 NT1030C-18M

